

ElectroCompetent Cells

Protocol:

- (1) Streak *E. coli* strain onto LB plate.
- (2) Pick O/N into 5 ml LB.
- (3) Inoculate 1.0 ml of the O/N into 1000ml LB.
- (4) Grow to mid log phase ($OD_{590} \sim 0.600$; about 3-5 hours). (To determine the cells' OD_{590} use the DU640 spec.)
- (5) Chill cells by placing the flask into an ice bath.
- (6) Spin cells down at 2,500 X g (~4200 RPM in SLA1500 rotor) for 5min at 4°C.
- (7) Resuspend cells in 1000ml ice cold sterile dH₂O.
- (8) Spin cells down 10,000 X g (~8500 RPM in SLA1500 rotor) for 15min at 4°C.
- (9) Repeat with a second 1000ml volume of ice-cold sterile dH₂O.
- (10) Repeat with 500ml ice-cold sterile dH₂O.
- (11) Resuspend in 20ml ice-cold sterile 10% Glycerol.
- (12) Transfers cells to an Oak Ridge tube.
- (13) Spin cells down at 10,000 X g (~9000 RPM in SS34 rotor) for 15min at 4°C.
- (14) Resuspend in 2ml ice-cold sterile 10% glycerol.
- (15) Aliquot into convenient volume (50 - 125µl).
- (16) *Quick-freeze tubes in liquid N₂.
- (17) Store cells at -80°C.

Media and Stuff:

- | | |
|--|-----------|
| (1) LB (1000 ml) | |
| Bacto tryptone | 10.0g |
| Yeast extract | 5.0g |
| NaCl | 2.5g |
| NaOH | 2 pellets |
| dH ₂ O | to 1000ml |
| (2) Sterile dH ₂ O | ~2.5L |
| (3) 10% glycerol | |
| (4) 0.5 ml Microfuge tubes | |
| (5) 250 ml centrifuge bottles | |
| (6) Oak Ridge Tubes | |
| (7) Sterile pipettes (pre-chilled) | |
| (8) Sterile pipette tips (pre-chilled) | |
| (9) SLA1500 rotor (pre-chilled) | |
| (10) SS34 Rotor (pre chilled) | |

Notes:

All media, solutions, centrifuge tubes, microfuge tubes and things that come into contact with the cells need to be sterile.

I keep the NaCl concentration low in my LB when making *Electrocompetent* cells.

Cells must be kept on ice the entire time cold after harvest (work in cold room or keep cells on ice).

I wash my cells 2 X with equal vol dH₂O whereas other protocols only do it once.

Cells don't pellet as well after being washed in H₂O; hence the increase in g forces; alternatively the length of the centrifuge run may be increased.

*Optional; i.e., cells can be placed directly into the -80°C freezer without freezing them first in N₂.