

**Oligonucleotide labeling**  
(Random Hexamer Primed and/or Directed primed)

Procedure:

- |                                       |          |
|---------------------------------------|----------|
| (1) DNA (GI fragment)                 | 50-100ng |
| (2) ddH <sub>2</sub> O                | to 50μL  |
| (3) Primers                           | 3.0μL*   |
| (4) Boil                              | 10min    |
| (5) Chill on ice                      |          |
| (6) Solution A                        | 2μL      |
| (7) HEPES                             | 2μL      |
| (8) dA                                | 1μL      |
| (9) dG                                | 1μL      |
| (10) dT                               | 1μL      |
| (11) BSA                              | 2μL      |
| (12) <sup>32</sup> P-dCTP             | 5μL      |
| (13) Klenow                           | 1μL      |
| (14) 37°C                             | ~1h      |
| (15) STOP                             | 100μL    |
| (15) Run G50 column (gravity or spin) |          |

Stock solutions and stuff:

(1) Primers (random hexamers) 1.8 mg/ml in T<sub>10</sub>E<sub>1</sub> (90 A<sub>260</sub> Units/ml)

(2) Sol'n A

- (a) Solution O 56μL
- (b) 2-ME 1μL

(3) Sol'n O

	<u>10ml</u>	<u>[Final]</u>
(a) Tris base	1.5g	1.25M
(b) MgCl	0.25g	0.125M
(c) pH	to 8.0	
(d) ddH <sub>2</sub> O	to 10ml	
(e) s/f		

(4) HEPES

	<u>10ml</u>	<u>[Final]</u>
(a) HEPES	4.8g	2M
(b) pH (NaOH)	to 6.6	
(c) ddH <sub>2</sub> O	to 10ml	
(d) s/f		

(5) dNTP 10mM stocks (make up according to Maniatis)

(6) BSA 10mg/ml (I use the stuff sent w/restriction enzymes)

(7) <sup>32</sup>P-dCTP (3000 Ci/mMole; NEN 013H)

(8) Klenow ~5U/μL

(9) Stop

	<u>10ml</u>	<u>[Final]</u>
(a) 5M NaCl	40μL	20mM
(b) 1M Tris (7.6)	200μL	20mM
(c) 0.5M EDTA	40μL	2mM
(d) 10% SDS	250μL	0.25%
(e) Orange G	3mg	0.03%
(f) Blue Dextran	100mg	1.0%
(g) ddH <sub>2</sub> O (sterile)	to 10ml	

Notes:

All stock solutions are kept in the -20°C.

Use either low melt agarose and β-agarase or a commercial kit to gel isolate template DNAs.

\*For directed priming reactions I use 3μL of a 125ng/3uL stock.

Make Solution A fresh every 3 months or so.

Ref:

Feinberg and Vogelstein (1983) A technique for radiolabeling DNA restriction endonuclease fragments to high specific activity. Anal Biochem. 132: 6-13.

Hodgson and Fisk (1987) Hybridization probe size control: optimized oligo labeling. Nuc. Acids Res. 15:6295